

# Tubigrip® Arthro-pad

## FOAM PADDED TUBULAR ELASTICATED SUPPORT BANDAGE

Tubigrip Arthro-pad provides non-invasive mechanical support for the elbow or knee after injury or surgery.



### BENEFITS OF TUBIGRIP ARTHRO-PAD

- Provides firm, comfortable support for the affected joint
- Retains body heat
- Easily applied and removed
- Washable for repeated use

 Tubigrip® Arthro-pad



## HOW TUBIGRIP® ARTHRO-PAD WORKS

Tubigrip Arthro-pads are constructed from polyurethane foam enclosed in a double layer of Tubigrip tubular bandage.

This double layer of Tubigrip provides firm support, encouraging use of the affected joint. The combination of Tubigrip and polyurethane foam also retains body heat.

Arthro-pads are quick and easy to apply, making them suitable for outpatient and home use.

They are washable for repeated use. The natural color of Tubigrip Arthro-pad is cosmetically acceptable to most patients.

## AREAS OF USE

Tubigrip Arthro-pads are used to provide support for the elbow or knee after injury or surgery. They may also be of value to patients with arthritis of these joints.

## MATERIAL PROPERTIES

- Polyurethane foam
- Elasticated support bandage consisting of:
  - Cotton
  - Elastodiene (contains natural rubber latex)
  - Polyamide

## WASHING INSTRUCTIONS



### Tubigrip Arthro-pad Ordering Information

Packaged sterile in single packs.

Product Code	Joint Circumference	Size
1585	Up to 12" (30cm)	Small
1586	Up to 16" (40cm)	Medium
1587	Up to 20" (50cm)	Large
1588	Up to 24" (60cm)	X-Large
1589	Up to 28" (70cm)	XX-Large

**Caution:** This product contains natural rubber latex which may cause allergic reactions, including anaphylactic responses.

